



Imaging's next chapter:

Illuminated with the intensity of natural light

> Rahel A. Kubik-Huch


Kantonspital Baden, Switzerland

rahel@ksb.ch



Clinical / Industry Relationship

- 8 years value partnership KSB SHS
- SOMATOM X.ceed Clinical Use Test
- Goals of translational projects → Combining strengths of both partners
 - Validation Studies: real-world patient data
 → improve SW & HW
 - Regulatory Clearance: studies as part of regulatory clearance (CE & FDA)
 - Market launch: supporting marketing activities, strengthening each others reputation





Translational studies SOMATOM X.ceed

- Interventional Radiology with an integrated 3D laser
- Advantages of wide bore, e.g., IR procedures, bariatric imaging
- Improved Diagnostic Image Quality, High-Res & Low-Dose








Assessment of Image Quality in Chest CT With Precision


Assessment of Image Quality in Chest CT With Precision
Matrix and Increased Framing Rate Using Single Source CT:
A Phantom Study


STYLIANOS MAVRIDIS¹, MONA EL-GEDAILY¹, RAHEL A. KUBIK-HUCH¹, FRIEDRICH KNOTH², JESUS FERNANDEZ LEON², ANDRÉ EULER³, LUKAS HEFERMEHL⁴ and TILO NIEMANN¹


> ¹Department for Radiology, Kantonsspital Baden, Baden, Switzerland; ²Siemens Healthcare GmbH, Forchheim, Germany; ³Department of Radiology, University Hospital Zurich, Zurich, Switzerland; ⁴Department of Urology, Kantonsspital Baden, Baden, Switzerland





doi:10.21873/invivo.13058













Gantry Integrated Guidance System for CT Interventions







Article

Evaluation of Accuracy and Performance of a Novel, Fully Gantry Integrated 3D Laser System for Computed Tomography Guided Needle Placement: A Phantom Study


Katharina Weigel ^{1,2}, Rahel A. Kubik-Huch ¹, Matthias Baer-Beck ³, Christian Canstein ³ and Michael Kostrzewa ^{1,4,*}











Gantry Integrated Guidance System for CT Interventions

- Interventions: faster, fewer needle corrections, less radiation dose
- Similar accuracy in needle placement novice & expert
- Full gantry integration avoids set up times → increases usage

Constant user feedback lead to system improvement:

- Improvement in HW & SW reliability
- Improvement in GUI
- Feedback towards development of dedicated HW & SW for IO




Figure 5. Procedural times for conventional needle guidance (CNG) and 3D laser guidance (3DLG), for expert and novice in in-plane and out-of-plane procedures.







Key Outcomes & Benefits

- Collaboration drives Innovation
 - Knowledge transfer: combining industry expertise with academic
 & clinical knowledge
 - Pooling of resources: access to specialized equipment, personnel & patients
 - Acceleration of progress: real-time feedback & improvements
- Collaboration benefits clinical institutions, industry & patients

Advancing precision medicine for patients' benefit





Thanks to

PD Dr Michael Kostrzewa, Prof. Tilo Niemann, Prof. André Euler & Team Kantonsspital Baden & the collaboration with SHS – Jesús Fernández León & CT Team